

# p!nto

# Mold Seat

for Wheelchair

Japanese ergonomic cushion for wheelchair




**A good seat gives you energy for an active life.**

P!NTO Mold Seat stimulates the potential of people in wheelchairs, it helps removing any tensions to make the use of your body easier. When you sit on P!NTO, feel how this cushion fits your body. The lateral supports maintain their mobility, the pelvis and the chest are supported. P!NTO reduces the mechanism of fatigue, eliminates stress, as a good seat gives you energy for an active life.

Designed in Japan, Production in Japan, Distribution from Italy.



## 4 ideal integrated supports



### Comparison of the zones of tension according to the posture



Just placing the P!NTO cushion on the wheelchair offers a comfortable and relaxing seat. You feel that your body is relaxed and your head is lighter. This cushion maximizes everyday comfort.

The body is constantly moving while you are sitting. The P!NTO cushion creates a position in which the body can return to its natural axis, relieving fatigue.

It adapts to different types of wheelchairs (with high back or adjustable back), can also be used on chairs, sofas, car seats and in every situations of the day.

This cushion prevents the body from bending forward and reduces pressure on the buttocks area. On a wheelchair without pinto, the instability of the posture causes tension in the shoulders and arms, limits comfortable movements and leads to tiring out. The pelvis tends to bend backwards and induces the user to unconsciously apply inverse tension to avoid slipping, possibly leading to stress and pain.

Thanks to this cushion, you will feel energized and up for the day.

### Technical characteristics:

Cushion Material: High Resilience HR Polyurethane Foam

Cover: 100%  $\mu$ -func® (mixed thread constituted by pure silver wool and special transparent polyester)

Dimensions, Weight: width 40cm x height 40cm x depth 40cm, 960g

Resistance: tested according to the JIS standard (Japanese Industrial Standard), thickness loss of only 4% after 80000 compressions

Maintenance: Cover can be removed and machine washed. The inner cushion does not require maintenance.

Portable: the polyurethane cushion, light and foldable can be easily transported

Use: thanks to the anti slip finish of the bottom side the cushion stays firm even when you move